## CONTROL UNIT FOR bioMIXdrive 1 STIRRING DRIVES

Designed for bioMIXdrive 1 stirring drive. For chemical, biotechnical and medical applications of all kinds.



## bioMIXcontrol 4MS

Efficient control unit, for 4 drives 2mag bioMIXdrive 1, developed for the hard daily lab use and uncompromising warming-free continuous operation.

- № Speed range 5 250 rpm
- ₯ For 4 stirring drives bioMIXdrive 1
- № Outstanding long-life quality
- Neveloped und Made in Germany







Туре	bioMIXcontrol 4MS
Speed range	5 - 250 rpm (4x individual settings)
Stirring power (max.)	6 W / bioMIXdrive 1 (automatic)
Power setting	-
Connection sockets	4
Material housing	stainless steel
Permitted operation conditions	0 °C up to +40 °C (at 80% humidity)
Electrical data	100-240 V / 50-60 Hz / 1.5 A
Dimensions (WxDxH)	200 x 155 x 38 mm
Protection category	IP20
Weight (gross)	approx. 2.0 kg
Order no.	98604

3 years warranty on material and manufacturing, reserve the right to make technical alterations. We do not assume liability for errors in printing or typeing.

Control unit for the connection of 4 stirring drive bioMIXdrive 1, individual and independent speed setting for each stirring drive (Master-Slave-Technology), also ON/OFF function, synchronous speed mode for all drives, extended speed range from 5 to 250 rpm, 6 watts/stirring position without any own-generated heating effect, automatic stirring power for the warming-free continuous operation e.g. in CO<sub>2</sub>-incubators, jerk-free stirring, even and smooth stirring also at slow speeds.

Digital display for speed setting, SoftStart for reliable and safe accelerating of stirring bar, QuickSet for quick setting of start and maximum speed, storage of last working status when turning off, firm stainless steel housing, easy to clean.

Compatible stirring drives bioMIXdrive 1



www.2mag.